

## **DNA Extraction** **(with Qiagen *DNeasy Plant Kit*)**

1. To dried mycelium, add:

400  $\mu$ l of AP1  
2  $\mu$ l of Rnase A
2. Grind solution with power drill and vortex to mix.
3. Incubate the mixture for 10 min at 65°C. Mix 2-3 time during incubation by inverting tube.
4. Add 130  $\mu$ l of Buffer P3 to the lysate, mix and incubate for 5 min on ice.
5. Centrifuge the lysate for 5 min at 13,000 rpm.
6. Transfer clear liquid to PURPLE QIAshredder Mini Spin Column placed in a 1.5 ml microcentrifuge tube. (The cap may need to be cut off so the tube will fit in the centrifuge.) Centrifuge for 2 min at 13,000 rpm.
7. Transfer the flow through fraction from the above step to a new tube and discard the Purple Spin Column. There should be about 400  $\mu$ l.
8. Add 1.5 the volumes of Buffer AW1 to the cleared lysate and mix by pipetting.

Ex. 400  $\mu$ l  $\rightarrow$  600  $\mu$ l of AW1 buffer
9. Pipet 500  $\mu$ l of the solution to the CLEAR Dneasy Mini Spin Column in a 1.5 ml microcentrifuge tube (with the cap cut off). Centrifuge for 1 min at 8000 rpm. Discard flow through.
10. Repeat Step 9 with remaining sample. Discard flow-through portion and collection tube.
11. Place the Clear Spin Column in a new 1.5 ml tube, add 500  $\mu$ l of Buffer AW2 to the Spin Column. Centrifuge for 1 min at 8000 rpm.
12. Discard flow-through, add 500  $\mu$ l of Buffer AW2 to Column and centrifuge for 2 min at 13,000 rpm. Centrifuge again for 2 min at 13,000 rpm to make sure there is no alcohol residual.
13. Transfer the Spin Column to a new 1.5 ml microcentrifuge tube (to be used to store the DNA, so label the cap), pipet 50  $\mu$ l of Buffer AE directly into the DNeasy membrane.
14. Incubate for 5 min at room temperature and then centrifuge for 2 min at 13,000 rpm.